



How to Setup your Reemo Watch

Setup Instructions v1.2

1 Make sure the watch is charged by placing it in the charging stand. Charging stand setup instructions are in the Quick Start Guide.



2 Press and hold the Power Button (on lower right) to turn on the watch.



3 Go to my.reemohealth.com/signup



To setup your account and activate your watch, you'll be asked for:

Account Information:

Name, Phone, Email, and a new Password for login.

Wearer Information:

Home Address, Date of Birth, and Gender.

Emergency Contacts:

Name, Email, Phone, and Relationship of up to two contacts. At least one contact is required, additional contacts can be added after registration.

If you need assistance, please visit reemohealth.com/support